

# Pandemic Flu Disasters Emergencies

**Be Prepared**  
with the Newport HT50 Ventilator



## The Newport HT50 Ventilator

### Clinically Effective

- Up to 10 hours of operation from internal battery
- Pediatric to Adult application
- Noninvasive or invasive
- Uses a simple, non-proprietary, inexpensive circuit for all patients
- Built-in gas generator eliminates need for external gas supply
- Built-in PEEP  
Pressure Support  
Volume Monitoring  
Back-up Ventilation  
Comprehensive Alarms

### For Emergency Preparedness

- Meets AARC Guidelines for "Ventilators to Meet Demands for Pandemic Flu and Mass Casualty Incidents"
- Rugged and lightweight
- Compatible for use with NATO CBRN filter to protect patients from chemical, biological, radiological and nuclear inhalation agents
- CEU programs available for Emergency Management Clinician training
- 24/7 Clinical and Technical Support



*Newport*  
**MEDICAL**

1620 Sunflower Ave., Costa Mesa, CA 92626  
800.451.3111 • www.ventilators.com