

Recharging the Dual Pac[™] Internal Battery System

Dual Pac[™] Internal Battery System - Longer Battery Use Time and Fast Recharge Time

- Primary lead acid battery provides up to 10-hours of use time depending on ventilator settings and battery condition.
- Secondary NiMH battery provides a minimum of 30 minutes use time after the Battery Low Alarm is violated.
- Both internal batteries simultaneously charge from any external AC or DC power.
- Automotive lighter DC power cord available.
- Fully depleted Dual Pac[™] Internal Battery System will re-charge to 100% in 7 to 8 hours.



Hospital Grade
Power Cord

Automotive Lighter Cable
(p/n HT460085)

**HT50 gives your patients the power to
Go and Go and Go!**