

**Setting up  
Newport HT50<sup>®</sup>  
Ventilator with  
a Reusable  
Heated Wire  
Single Limb  
Circuit**

The following items are needed in order to set up the Newport HT50 Ventilator with a reusable heated wire, single limb circuit for use with a Fisher & Paykel humidifier.

**→ For use with the Fisher & Paykel MR850 Humidifier**

(Note that this humidifier requires use of a heater wire in the circuit.)

P/N	Description
HT600045	Reusable breathing circuit with exhalation valve
900MR858	Heater wire adapter for reusable circuits
900MR510	Heater wire assembly for 52" (1.3 meter) circuit tubing
900MR860	Temperature probe for 52" (1.3 meter) circuit tubing
TUB024P	Circuit tubing for connecting the HT50 <sup>®</sup> with a 3 <sup>rd</sup> party humidifier
900MR070	60" (1.5 meter) Draw wire

**→ For use with the Fisher & Paykel HC500\* Homecare Humidifier**

(Note that this humidifier can be used with or without a heater wire in the circuit.)

P/N	Description
HT600045	Reusable breathing circuit with exhalation valve
900MR558	Heater wire adapter for reusable circuits
900MR510	Heater wire assembly for 52" (1.3 meter) circuit tubing
900MR560	Temperature probe for 52" (1.3 meter) circuit tubing
TUB024P	Circuit tubing for connecting the HT50 <sup>®</sup> with a 3 <sup>rd</sup> party humidifier
900MR070	60" (1.5 meter) Draw wire

\* Order Directly from Fisher & Paykel.