


TO: Newport Customers

SUBJECT: Replacement Parts and Accessories Available for the Expiratory Filter Heater Kits KIT2110A & KIT2220A

The purpose of this bulletin is to inform you of the part numbers for ordering replacement parts and accessories for the Expiratory Filter Heater Kits KIT2110A & KIT2220A.



Tubing:
 Reusable 22 mm ID Hytrel 18" (45 cm) Circuit Tubing (p/n TUB018P)
 Reusable 22 mm ID Hytrel 24" (60 cm) Circuit Tubing (p/n TUB024P)
 Reusable 22 mm ID Silicone 18" (45 cm) Circuit Tubing w/Elbow (p/n 9013-01-1)

Filter:
 Reusable Filter (p/n FLT1200P)
 Disposable Filter (p/n 3000/03-C)

Filter Heater Ring (F-4030-01)

Heater (does not include mounting hardware)
 220-240 v Heater (p/n 800-22202)
 100-120 v Heater (p/n 800-21101)

Silicone Connector (p/n ADP1802P)

Reusable 75 ml Water Trap (p/n G-313003-1)

Additional Heater Accessories:

- **Dovetail (male) Mount Bracket (p/n BKT1808M):**
REQUIRED for mounting the heater to the cart, the Rail Mount Bracket (p/n 900MR088) or Dual Rail Mount Bar (BAR1820A).
- **Dual Rail Mount Bar (p/n BAR1820A):**
REQUIRED for mounting both a Filter Heater and an F&P humidifier to cart.
- **Rail (Female) Mount Bracket (p/n 900MR088):**
 For an easy quick-release mounting to standard e500 cart or e360 Fairfield Rail.



Additional Heater Accessories Cont'd:

- **Disposable Water Trap Kit (p/n WTR134P):**

Includes suction tubing for attachment to a suction source for safe removal of condensate.



Bacteria Filter Heater Parts and Accessories:

P/N	Description
KIT2110A	Bacteria Filter Heater Kit 110V w/ reusable: water trap, filter, circuit tubing and silicone connector (requires BKT1808M)
KIT2220A	Bacteria Filter Heater Kit 220V w/ reusable: water trap, filter, circuit tubing and silicone connector (requires BKT1808M)
TUB018P	Reusable 22 mm ID Hytrel 18" (45 cm) Circuit Tubing
TUB024P	Reusable 22 mm ID Hytrel 24" (60 cm) Circuit Tubing
9013-01-1	Reusable 22 mm ID Silicone 18" (45 cm) Circuit Tubing w/Elbow
FLT1200P	Reusable Bacteria Filter
3000/03-C	Disposable Bacteria Filter (12/case)
F-4030-01	Filter Heater Ring
800-21101	Filter Heater Only 100-120 V
800-22202	Filter Heater Only 220-240 V
WTR134P	Disposable Water Trap
G-313003-1	Reusable Water Trap (75 mL)
ADP1802P	Reusable 22 mm ID Silicone Connector
TUB018P	Reusable 22 mm ID Circuit Tubing (18" / 45.7 cm)
ADP1802P	Reusable 22 mm ID Silicone Connector, to connect water trap to heater
BKT1808M	Male Dovetail for back of Filter Heater to allow it to be mounted to cart, dual bar, etc.
900MR088	Rail Mount Bracket to Mount Filter Heater onto e500 Cart or Fairfield rail
BAR1820A	Dual mounting bar with two female dovetail mounts on it. To mount both a humidifier and a Filter Heater onto a Newport e360 or e500 Cart.

To order, please contact Newport Customer Service at ++1.714.427.5811 ext. 282 or email: customers@ventilators.com.

MKPB042110