


# PRODUCT BULLETIN

24 February, 2003

**TO: NEWPORT Distributors**  
**SUBJECT: Fisher-Paykel Mounting Tray for use with Newport HT50**



Some customers may be looking for a way to secure equipment at the bedside. The Mounting Tray accommodates the Newport HT50, F&P humidifier and Water Bag Pole on a stable platform. This helps prevent accidental spills or circuit disconnects and helps organize the patients bedside equipment.

The Water Bag Pole is optional and is ordered separately. It screws into the top of the dove tail mount for the humidifier.

**New HT50 Accessory**  
**Mounting Tray**  
**Water Bag Pole**

**Part Number**  
**900HC032**  
**900MR012**



The picture at left illustrates the complete setup on a night stand. The Mounting Tray is 20" by 12" and weighs just over 2 pounds. The Water Bag Pole is 3' 6" tall. The Mounting Tray and Pole together weigh just under 3 pounds.

The Mounting Tray comes with a Velcro strap that secures the HT50 to the tray. If the HT50 needs to be moved to a wheel chair for transport, it can be easily disconnected from the tray.



If you have any questions or need additional information please contact NEWPORT Customer Service: 1.949.642.3910 ext 282 or email: [info@newportnmi.com](mailto:info@newportnmi.com).

MKPB022403

**NEWPORT  
ENMIE**